



Non-Woven Roll with Pad (Instructions for Use)

The Asguard Non-Woven Roll with Pad combines an absorbent pad with a soft and comfortable non-woven adhesive backing. The self-adhesive rolls are "cut-to-length" and are ideal for fixation of difficult contour areas. The adhesive is skin-friendly and provides reliable fixation, yet is easy to remove, minimising trauma to healing tissue.




STEP 1

Cleanse wound thoroughly, using local protocols. The skin surrounding the wound should be clean and dry.




STEP 2

Dressing size is dependent on the wound surface area. Use grid to determine length – the dressing edge should extend 2-3cms beyond the wound margin.




STEP 3

Cut to desired length. Fold both backing papers to expose dressing pad.




STEP 4

Apply the Asguard Non-Woven dressing covering entire wound bed. Release backing paper.



STEP 5

Smooth down the edges of the non-woven border ensuring adherence.



DRESSING CHANGE

The clinician-in-charge is responsible for determining the need for dressing changes, dependent on the stage and phase of wound healing and exudate level.

The dressing should be changed when exudate absorption is considered to be extending beyond the wound surface area.

7 days is the maximum period between dressing changes.

To minimise potential damage to fragile skin, Asguard Non-Woven Roll with Pad should be eased off using an adhesive remover wipe.

INDICATIONS

Asguard Non-Woven Roll with Pads are non-sterile and indicated for use only as primary dressing in a first-aid or first response situation in the management of:

- Abrasions and lacerations
- Blisters and minor skin loss

CONTRAINDICATIONS

Should not be used if allergies to any of its components, is known.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

Sentry Medical Pty Ltd

1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au