



ISP[®] Dual Steam Monitoring Strip



ISP[®] Dual Steam Monitoring Strip

The ISP[®] Dual Steam Monitoring Strip is a multi-variable indicator designed to indicate whether or not the lower limits of the essential process parameters have been achieved throughout the sterilization cycle.

The ISP[®] Dual Steam Monitoring Strip has a perforation in the middle. When separated, one part can be stored in the patient file and the other part of the strip can be returned to the CSSD to complete the sterilization file.

The non-toxic indicator ink will change from light blue to an even dark colour under the correct sterilization conditions, i.e. steam penetration time and temperature. The strips are laminated to prevent staining of the ink on the instruments or products.

The ISP[®] Dual Steam Monitoring Strip complies with EN-ISO 11140-1 (Class 4).



TEL +612 9609 5995
FAX +612 9609 5991
EMAIL sales@sentrymedical.com.au
WEB www.sentrymedical.com.au

Art. No.: 33002014
Quantity: 250 pieces/box

INTERSTER [®]
Interster International BV