


ISP[®] Helix Test System


ISP[®] Helix Test System

The Helix Test System by Interster is the ideal process challenge device to perform a “worst case” scenario when sterilizing instruments. This test is suitable for every sterilization process in a loaded autoclave. It can also be used to ensure maximum steam penetration where hollow instruments, or instruments with small lumina, are concerned.

The Helix strips are in compliance with EN 867-5 and EN-ISO 11140-1 (Class 6). It is easy to interpret the strip due to the clear and even colour change. The ink used is non-toxic and it can therefore not cause any damage to instruments. Due to the slice strip at the back side, the strip can be stored easily in a log or on a registration card.



How to use the test:



Art. No.: 33002030
Quantity: 100 pieces/box



TEL +612 9609 5995
FAX +612 9609 5991
EMAIL sales@sentrymedical.com.au
WEB www.sentrymedical.com.au

INTERSTER[®]
Interster International BV