



Low
Linting



Available
in 3 Grades



Dual
Coloured

Sterifit Sterilisation Wraps are made from SMS bonded polypropylene. They are suitable for steam, ethylene oxide, and plasma sterilisation processes. Dual colouring allows for easy detection of strike-through, tears and punctures.

DUAL COLOUR (BLUE/GREEN)


For fast and efficient detection of tears

STRONG AND DURABLE

Provides a high resistance to tearing and abrasion

SUITABLE FOR ALL STERILISATION METHODS

Can be used in steam, ethylene oxide, or plasma sterilisation processes



Suitable for:

STANDARD

- Light to medium tray
- Linen packs
- Medium devices & instruments
- Basin sets


HEAVY DUTY

- Large trays
- Instrument sets
- Medical devices
- Multiple basin sets

ULTRA HEAVY DUTY

- Heavy surgical sets such as orthopaedic
- Heavy instrument sets

*Reorder details on back



Reorder	Description	Size	Weight	Colour	Sterility	Pcs/Unit	Units/Ctn
STANDARD							
STS102A	Sterifit 2 - Standard (Dual)	60cm x 60cm	104gsm	Green/Blue	Non-Sterile	250	1
STS103A	Sterifit 2 - Standard (Dual)	76cm x 76cm	104gsm	Green/Blue	Non-Sterile	150	1
STS104A	Sterifit 2 - Standard (Dual)	91cm x 91cm	104gsm	Green/Blue	Non-Sterile	75	1
STS109A<<	Sterifit 2 - Standard (Dual)	101cm x 101cm	104gsm	Green/Blue	Non-Sterile	75	1
STS106A	Sterifit 2 - Standard (Dual)	121cm x 121cm	104gsm	Green/Blue	Non-Sterile	50	1
HEAVY DUTY							
STH101A	Sterifit 2 - Heavy Duty (Dual)	91cm x 91cm	122gsm	Green/Blue	Non-Sterile	75	1
STH102A	Sterifit 2 - Heavy Duty (Dual)	101cm x 101cm	122gsm	Green/Blue	Non-Sterile	75	1
STH103A	Sterifit 2 - Heavy Duty (Dual)	121cm x 121cm	122gsm	Green/Blue	Non-Sterile	25	1
ULTRA HEAVY DUTY							
STU102A	Sterifit 2 - Ultra Heavy Duty (Dual)	121cm x 121cm	132gsm	Green/Blue	Non-Sterile	25	1
STU103A<<	Sterifit 2 - Ultra Heavy Duty (Dual)	137cm x 183cm	132gsm	Green/Blue	Non-Sterile	25	1

<< Indent Item