



DUST COVER BAGS



Constructed from low density polythene, making them highly durable. These bags can be heat sealed. Dust Cover Bags provide an ideal storage environment for instruments and trays that have undergone the sterilisation process, and need protection from dust whilst in storage.

PRODUCT FEATURES

- o Durable
- o Can be heat sealed

SUITABLE FOR

- o Providing storage for sterilised equipment

REORDER	SIZE	COLOUR	PCS/UNIT	UNITS/CTN
DC001	15.2cm x 25.4cm	Clear	100	20
DC002	20.3cm x 30.5cm	Clear	100	15
DC003	25.4cm x 38.1cm	Clear	100	15
DC004	30.5cm x 50.8cm	Clear	100	15
DC005	40.6cm x 50.8cm	Clear	100	5
DC006	45cm x 76.2cm	Clear	100	4
DC007	61cm x 76.2cm	Clear	100	4

All product information is correct at date of publication

