

SELF SEALABLE DUST COVER BAGS



The ISP® Self Sealable Dust Cover Bags are constructed from a polyester and polypropylene blend. Dust Cover Bags provide an ideal storage environment for instruments and trays that have undergone the sterilisation process, and need protection from dust whilst in storage.

PRODUCT FEATURES

- o Easy peel opening
- o Adhesive sealing tab with area for clear labelling
- o Conform EN-ISO 11607-1

SUITABLE FOR

- o Storage of sterilised instruments and trays
- o Extends sterile tenability

REORDER	SIZE	COLOUR	PCS/UNIT	UNITS/CTN
3FSSV510802	15cm x 25cm	Clear	100	5
3FSSV510804	20cm x 30cm	Clear	100	5
3FSSV510806	25cm x 37.5cm	Clear	100	5
3FSSV510808	30cm x 50cm	Clear	100	5
3FSSV510810	40cm x 55cm	Clear	100	5
3FSSV510812	40cm x 75cm	Clear	50	5
3FSSV510814	60cm x 75cm	Clear	50	5