



All Burn
Types



Suitable for
Children



Ideal for
First Aid

Burnshield is a sterile trauma hydrogel which can be used for minor burns and scalds. With its high water content, Burnshield Hydrogel cools and soothes by absorbing and dissipating heat which minimises skin damage. The gel is clear which makes it easy to assess the wound. It is non-stick and will not adhere to the wound. It is non-toxic and non-irritant. Safe to use on children.

SAFE FOR USE ON CHILDREN

Non-toxic and non-irritant formula is safe for all ages

ABSORBS AND DISSIPATES HEAT

Gel structure quickly cools and soothes the affected area minimising burn damage



Suitable for:

Superficial / partial thickness burns

Rehydration of dressings

Difficult to reach areas

All types of burns: fire, sunburn, thermal & scalding, chemical, radiation

Reorder	Description	Size	Pcs/Unit	Units/Ctn
550004	BURNSHIELD® Burn Sachets	3.5ml	1200	1
550009	BURNSHIELD® Hydrogel Tube	25ml	72	1
550000	BURNSHIELD® Hydrogel Bottle	50ml	72	1
550020	BURNSHIELD® Hydrogel Spray	50ml	48	1
550003	BURNSHIELD® Hydrogel Spray	125ml	48	1
550065<<	BURNSHIELD® Blott Blister Pack (5pk)	5 x 3.5ml	100	1
550016<<	BURNSHIELD® Hydrogel Tube & Dispenser	12 x 25ml	12	1

<< Indent Item. MOQ required