



Protects
Instruments



Various
Sizes



Easy
To Use

Interster Instrument Protectors are ideal to protect sharp and fragile instruments, the sterilisation packaging material and to protect the user. A batch number is printed on each protector for full traceability of critical surgical instruments.

TRACEABLE

Batch number printed for full traceability

CLEAR FILM POCKET

Protects sharp tips and edges

MULTIPURPOSE

Steam, ethylene oxide, formaldehyde sterilisation compatible



Suitable for:

Small critical surgical instruments

Optimal protection for instruments

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
3FSZB730802	Small	Non-Sterile	100	1
3FSZB730804	Medium	Non-Sterile	100	1