



Leak-Proof



Easy  
To Use



Easy  
Inspection

Interster Liquid Control Dye Penetration Tests are specially developed to check the integrity of the seals, produced with different heat sealers. You can safely do your own quality dye test without damaging the seal. Other advantages are that you can prepare and seal the test packaging in advance, without the risk of leakage or migration of blue or red liquid on valuable objects in your work area.

## NO PIPETTES USED

Safely perform test without damaging the seal

## PREPARE IN ADVANCE

Removes risk of leakage



## Suitable for:

Type 1 - See-through packaging which consist of a film/medical grade paper combination

Type 2 - See-through packaging which consist of a film/non-woven polyolefins (Tyvek®) combination

Reorder	Type	Feature	Ink	Pcs/Unit	Units/Ctn
3FSZB830816	Type 1	Film / Medical Paper	Blue	30	1
3FSZB830818	Type 2	Film / Tyvek	Red	30	1