



Easy To Use



Clear Seal Indication

The Interster Seal Checks display a suitable seal and show indications of possible sealing problems for sterilisation pouches.

PROVIDES VALIDATION TO ENSURE

- Intact seal for a specific width
- No channels or open seals
- No punctures or tears in seals
- No material delamination or separation



Suitable for:

Paper/film and/or Tyvek sterilisation pouches or reels in combination with any rotary heat sealer

Reorder	Description	Feature	Sterility	Pcs/Unit	Units/Ctn
3FSZB830802	Seal Check	Paper / Film	Non-Sterile	250	1
3FSZB830803	Tyvek Seal Check	Tyvek / Film	Non-Sterile	125	1