



Writable Surface



Type 1 Process Indicator



Secure Adherence

The Interster Steam Indicator Tapes are intended for sealing packages to be steam sterilised. They indicate exposure to saturated steam as a Type 1 indicator. The non-toxic ink indicator has a light to dark colour change.

COLOUR CHANGE INDICATOR

Verifying exposure to sterilisation process

LEAD FREE

No hazardous waste

CAN BE WRITTEN ON

Surface can be written on without smudging



Suitable for:

Use on linen/woven material

Sealing packages to be steam sterilised

Reorder	Size	Colour	Sterility	Pcs/Unit	Units/Ctn
3FKLB410802	19mm x 50m	Buff	Non-Sterile	1	48
3FKLB410806	19mm x 50m	Blue	Non-Sterile	1	48
3FKLB410804	25mm x 50m	Buff	Non-Sterile	1	36
3FKLB410808	25mm x 50m	Blue	Non-Sterile	1	36