



Various Colours



Dries Instantly



Non-Toxic

Interster Sterilisation Markers are designed for light touch writing to avoid damage to packaging material. The solid writing point of the marker ensures a constant writing quality and long durability. The ink dries immediately and does not fade. The pen is for marking sterilisation packaging material, with special autoclave low odour ink, without additor or oxylene/toluene. The ink is steam, gamma, formaldehyde, plasma and EO resistant.

NON-TOXIC INK

No damage to instruments

INSTANTANEOUS DRYING INK

Neutral odour, dries immediately and does not fade

MULTIPURPOSE

Ink is steam, gamma, formaldehyde, plasma and EO resistant

Suitable for:

For marking sterilisation packaging material

Reorder	Colour	Sterility	Pcs/Unit	Units/Ctn
3FSZB710824	Black	Non-Sterile	25	1
3FSZB710826<<	Blue	Non-Sterile	25	1
3FSZB710828	Red	Non-Sterile	25	1

<< Indent Item

