



# Absorbent Medical Towel

AAXIS® Absorbent Medical Towel is a general purpose towel designed to be used as a surgeons scrub towel, hand towel or tray liner. Due to their versatility, AAXIS® Absorbent Medical Towels can also be used as a general wiper or bowl spacer.



## Features

Low linting

Soft and gentle on skin

Smart dispenser design and easy tear off

Absorbs 7 times its weight

Wet and dry tear resistance

Minimises risk of cross contamination

Latex-free

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10030400	AAXIS® Absorbent Medical Towel 30x35cm, 100's	12	N/A	12
10030401	AAXIS® Absorbent Medical Towel 30x50cm, 100's	12	N/A	12