



CPE Impervious Splash Gown

Features

- Impervious polyethylene material
- Universal size thumb loop
- Perforated upper back seams to allow easy doffing
- Waterproof heat welded arm seams
- Open back to reduce heat and improve breathability
- Lint-free, latex-free
- Non-sterile, single use
- Convenient inner box of 10 gowns



Waterproof seams

Universal thumb loop



Perforations for single handed doffing



PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/Inner	SKU/CTN
10025517	Aaxis® Splash Gown, Impervious CPE	100	10	100