




# AAXIS® Surgical Mask AS Level 2 & Type IIR



## Features

High quality, 3-ply Surgical mask

Complies with:

AS 4381: 2015

EN 14683: 2019 Type IIR (120mmHg)

Convenient dispenser box of 50

Non-sterile

Test Requirement	EN 14683:2019 Type IIR	AS 4381: 2015 Level 2	Aaxis Surgical Mask
BFE %	≥98	≥98	Pass
Synthetic blood resistance	≥ 16.0 kPa >120 mmHg	>120mmHg	Pass
Differential pressure	<60 Pa/cm2 <6.0 mmH2O/cm2	<5.0 mmH2O/cm2	Pass
Microbial Cleanliness	≤30 cfu/g	n/a	Pass

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10025518	Aaxis® Surgical Mask AS Level 2 & Type IIR Tie 50's	30	N/A	30
10025519	Aaxis® Surgical Mask AS Level 2 & Type IIR Loop 50's	30	N/A	30