

Emesis Bag / Retail Shelf Pack

Emesis Bags have a firm twist and lock rim to provide easy control, and help prevent spillage after use. They isolate splash back and have a graduated indicative volume system to 1500ml, allowing measurement of contents in the bag.

On completion of use, the bag is easily disposed of with no mess. Simply twist the bag and lock into the keyhole slot on the rim.



Fuss Free Use

Thoughtfully designed for easy inclusion in first aid kits and medicine cabinets



Infection Control

Easy non-leak twist and lock seal for a no mess disposal



Convenient Dispenser

Handy wall dispenser available for hospital setting, or 3 pack dispenser for retail environments



- Suitable for:**
- First aid
 - Medical settings
 - General home use

REF	Description	Sterility	Pcs/Unit	Units/Ctn
VB002	Emesis Bag 1500ml (Red Ring)	Non-Sterile	50	10
10110723	Emesis Bag 1500ml (Yellow Ring)*	Non-Sterile	50	12
VB011	Emesis Bag Retail Shelf Pack (3pcs/pk - 6pks)	Non-Sterile	10	1
VBD001	Acrylic Dispenser (Suitable for VB002 and 10110723)	Non-Sterile	1	12

*10110723 Does not come in a dispenser box

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

