

Emesis Bag / Retail Shelf Pack

Emesis Bags have a firm twist and lock rim to provide easy control, and help prevent spillage after use. They isolate splash back and have a graduated indicative volume system to 1500ml, allowing measurement of contents in the bag. On completion of use, the bag is easily disposed of with no mess. Simply twist the bag and lock into the keyhole slot on the rim.



Fuss Free Use
Thoughtfully designed for easy inclusion in first aid kits and medicine cabinets



Infection Control
Easy non-leak twist and lock seal for a no mess disposal



Convenient Dispenser
Handy wall dispenser available for hospital setting, or 3 pack dispenser for retail environments



- Suitable for:**
- First aid
 - Medical settings
 - General home use

REF	Description	Size	Type	Sterility	Pcs/Unit	Units/Ctn
VB002	Emesis Bag	1500ml	Red Ring	Non-Sterile	50	10
10110723*	Emesis Bag	1500ml	Yellow Ring	Non-Sterile	50	12
VB011	Emesis Bag Retail Shelf Pack	1500ml	3pcs/pk - 6pks	Non-Sterile	10	1
VBD001	Acrylic Dispenser	140 x 140 x 350mm	Suitable for VB002 and 10110723	Non-Sterile	1	12

*10110723 Does not come in a dispenser box

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

VOMIT BAG

GENERAL WASTE

POLYBAG

CHECK LOCALLY

BOX

RECYCLE ME