

Purifil™ Prefilled Water Syringe

The Purifil Prefilled Water Syringe is specifically designed to fill latex and silicon Foley indwelling Catheter balloons in a single step. The syringe's luer slip is inserted into the catheter balloon port directly to enable filling. The syringe is graduated allowing the clinician to measure the amount of sterile water that is to be filled from the syringe.




Easy to Use

Ready for use, no preparation required



Single Use

Prepackaged in individual vacuum peel pouches



Easy-to-Read Graduations

Has conventional and easy-to-read graduations



Suitable for: - Filling latex and silicon Foley indwelling catheter balloons

REF	Size	Sterility	Pcs/Unit	Units/Ctn
10045000	10 mL	Sterile	100	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

