





AAXIS® Patient Modesty Gowns

The AAXIS® Patient Modesty Gowns are a uniquely designed patient gown that is comfortable, breathable and opaque. The gown is constructed using a high-quality SMS material to ensure patients feel as comfortable as possible during examination. The gown is constructed using a new automated production assembly process to improve consistent quality and the gown is equally suitable for a wide variety of patient dress applications such as routine examination, radiological scans and other surgical and medical applications where access under clothing is required.

Uniquely, the gown can be fitted either as a jacket or can be walked into with the ties part of the construction of the gown itself. The ties are easily detached to a point under the arms to secure the gown around the body. Wearing the gown as a jacket can allow more easy access for breast radiological examinations, without full removal.



High quality SMS material: opaque, comfortable and breathable



All seams are sonic welded ensure complete patient comfort



Ties are easily detached from either side of the gown before or after donning




Features

- Comfortable, breathable and opaque SMS material
- Generous size cuts
- Can be used with the opening at the front or the back
- Convenient tear away ties
- Single patient use
- Non-sterile

Applications


- Patient preparation for radiological scans
- Routine examinations
- Other hospital surgical and medical applications

Fitment/Donning




1

Walk into the gown with the opening at the back. Detach the ties from the bottom up.



2

Cross the ties across your back.



3

Knot the ties together across your front or back.

PRODUCT CODE	PRODUCT DESCRIPTION	SIZE	TOTAL LENGTH (mm)	SHOULDER WIDTH (mm)	OPAQUE	SKU/TU (MOQ)	SKU/Inner	SKU/CTN
10026000	Aaxis® Patient Modesty Gown, M	M	1060	720	Y	100	25	100
10026001	Aaxis® Patient Modesty Gown, X-Large	L/XL	1310	830	Y	100	25	100
10026002	Aaxis® Patient Modesty Gown, XX-Large	XXL	1500	900	Y	100	25	100