

Asguard[®]

ACUTE DRESSING RANGE





Asguard[®]

For over 30 years, Sentry Medical has been helping health care professionals provide exceptional care to their patients through a modern and quality driven portfolio.

A full range of dressings and tapes, ideal for low exuding wounds and difficult to dress areas.

Asguard Dressings have been specifically designed to provide a secure barrier to protect against potential infection and contaminants.

The extensive range comprises of various unique shapes, aiding to dress difficult areas (such as joints) and are ideal for minor wounds and first aid.

Asguard®

Transparent Film Dressing / Film Roll (Clear)

A transparent adhesive polyurethane film with high moisture vapour permeability which adheres to the skin and forms a flexible and comfortable layer. The dressing and rolls conforms to the skin and allows easy visibility to wound healing.



Asguard®

Transparent Film Island Dressing (Clear+)


Asguard Transparent Film Island Dressings and Film Rolls provides a low adherent and breathable absorbent pad permeable to moisture vapour. The waterproof polyurethane film assists with unnecessary dressing changes and helps protect from external contaminants.



Asguard®

Non-Woven Roll (Flex)

Asguard Non-Woven Roll is a "cut to length" self-adhesive, non-woven fabric fixation tape. Highly comfortable, this tape is an ideal fixation of primary dressings for contoured areas and can flex with everyday movement.



Asguard®

Non-Woven Island Dressing / Roll with Pad (Flex+)


Asguard Non-Woven Island Dressings and Roll with Pads combine an absorbent pad with a soft and comfortable non-woven adhesive backing. The self-adhesive rolls are "cut-to-length" and ideal for fixation of difficult contour areas.



Asguard®

Reinforced Skin Closures (SurgiStrip)

Asguard Reinforced Skin Closures allow earlier post-operative suture or staple removal to assist natural healing and minimise scarring. They are a painless alternative for the closure of minor incisions or skin lacerations, enabling fast and easy application, with no need for anaesthetic and eliminating needle puncture marks.





Asguard Dressing Features & Benefits:

Asguard®		CLEAR DRESSING	CLEAR ROLL	CLEAR+ DRESSING	FLEX ROLL	FLEX+ DRESSING	FLEX+ ROLL	SURGI-STRIP
WATERPROOF		•	•	•				
TRANSPARENT		•	•	•				•
BREATHABLE		•	•	•	•	•	•	•
CAN BE CUT TO SIZE			•		•		•	
ABSORBENT PAD				•		•	•	
NON-WOVEN					•	•	•	•
DATE LABEL		•		•				
CONFORMABLE		•	•	•	•	•	•	•

Asguard Dressing User Guide

Asguard®	CLEAR	CLEAR	CLEAR+	FLEX	FLEX+	FLEX+	SURGISTRIP
	Transparent Film Dressing	Transparent Film Roll	Transparent Film Island Dressing	Non-Woven Roll	Non-Woven Island Dressing	Non-Woven Roll with Pad	Reinforced Skin Closures
PARTIAL THICKNESS WOUND	✓		✓		✓		
LITTLE TO NO EXUDATE	✓		✓		✓		✓
MINIMAL TO LOW EXUDATE			✓		✓		✓
STAGE 1 PRESSURE ULCER	✓		✓				
STAGE 2 PRESSURE ULCER	✓		✓				
REINFORCE IV SITES & LINES	✓	✓		✓		✓	
DONOR SITES	✓		✓		✓		
MILD LACERATIONS	✓		✓		✓		✓
ABRASIONS	✓		✓		✓		
SECOND DEGREE BURNS			✓		✓		
AFFIX / COVER PRIMARY DRESSING	✓	✓		✓		✓	✓
WATERPROOF BARRIER	✓	✓	✓				
NECROTIC WOUNDS	✓		✓				
PRIMARY DRESSING	✓		✓		✓		✓
SECONDARY DRESSING	✓	✓		✓		✓	
SECONDARY PADDING LAYER				✓		✓	

REF	DESCRIPTION	SIZE	PCS/UNIT	UNITS/CTN
ASGUARD TRANSPARENT FILM DRESSING / FILM ROLL (CLEAR)				
ISB001	Asguard Transparent Film Dressing	4cm x 5cm	50	10
ISB002	Asguard Transparent Film Dressing	6cm x 7cm	50	10
ISB003	Asguard Transparent Film Dressing	10cm x 12cm	25	10
ISB004	Asguard Transparent Film Dressing	15cm x 20cm	25	10
ISB1025	Asguard Transparent Film Dressing	10cm x 25cm	25	10
ISB2030	Asguard Transparent Film Dressing	15cm x 20cm	25	10
ISB010	Asguard Transparent Film Roll	5cm x 10m	1	24
ISB011	Asguard Transparent Film Roll	10cm x 10m	1	12
ASGUARD TRANSPARENT FILM ISLAND DRESSING (CLEAR+)				
ISC001	Asguard Transparent Film Island Dressing	4cm x 5cm	50	10
ISC002	Asguard Transparent Film Island Dressing	6cm x 7cm	50	10
ISC0610	Asguard Transparent Film Island Dressing	6cm x 10cm	50	4
ISC0910	Asguard Transparent Film Island Dressing	9cm x 10cm	25	10
ISC0915	Asguard Transparent Film Island Dressing	9cm x 15cm	25	10
ISC003	Asguard Transparent Film Island Dressing	10cm x 12cm	25	10
ISC028	Asguard Transparent Film Island Dressing	10cm x 20cm	25	10
ISC025	Asguard Transparent Film Island Dressing	10cm x 25cm	25	10
ISC026	Asguard Transparent Film Island Dressing	10cm x 30cm	25	10
ISC0935	Asguard Transparent Film Island Dressing	10cm x 35cm	30	8
ISC004	Asguard Transparent Film Island Dressing	15cm x 20cm	25	10
ASGUARD NON-WOVEN ROLL (FLEX)				
ISD020	Asguard Non-Woven Roll	2.5cm x 10m	1	30
ISD021	Asguard Non-Woven Roll	5cm x 10m	1	24
ISD022	Asguard Non-Woven Roll	7.5cm x 10m	1	24
ISD023	Asguard Non-Woven Roll	10cm x 10m	1	12
ISD024	Asguard Non-Woven Roll	15cm x 10m	1	12
ISD025	Asguard Non-Woven Roll	20cm x 10m	1	6
ASGUARD NON-WOVEN ISLAND DRESSING / NON-WOVEN ROLL WITH PAD (FLEX+)				
ISD001	Asguard Non-Woven Island Dressing	7.5cm x 5cm	50	10
ISD002	Asguard Non-Woven Island Dressing	8cm x 10cm	50	10
ISD030	Asguard Non-Woven Island Dressing	8.25cm x 6cm	50	10
ISD003	Asguard Non-Woven Island Dressing	10cm x 15cm	20	10
ISD004	Asguard Non-Woven Island Dressing	10cm x 20cm	20	10
ISD006	Asguard Non-Woven Island Dressing	10cm x 25cm	25	10
ISD008	Asguard Non-Woven Island Dressing	10cm x 30cm	25	10
ISD011	Asguard Non-Woven Roll with Pad	5cm x 10m	1	24
ISD012	Asguard Non-Woven Roll with Pad	7.5cm x 10m	1	24
ISD013	Asguard Non-Woven Roll with Pad	10cm x 10m	1	12
ASGUARD HAEMOPLAST DRESSING (FLEX+)				
ISD101	Asguard Haemoplast Dressing	3.6cm Ø	100	20
ISD100	Asguard Haemoplast Dressing	4cm x 5.5cm	100	20
ASGUARD REINFORCED SKIN CLOSURES (SURGISTRIP)				
C-01	Asguard Reinforced Skin Closures	3mm x 75mm	50	20
C-02	Asguard Reinforced Skin Closures	6mm x 75mm	50	20
C-03	Asguard Reinforced Skin Closures	6mm x 38mm	50	20
C-04	Asguard Reinforced Skin Closures	6mm x 100mm	50	20
C-05	Asguard Reinforced Skin Closures	13mm x 100mm	50	20
C-06	Asguard Reinforced Skin Closures	25mm x 125mm	25	20