

IV Dressing Reinforced

The Asguard Reinforced IV Dressings are constructed from a High MVTR film, which allows for excess moisture to evaporate, combined with a laminated non-woven border for securement and enhanced fixation. They maintain firm adhesion and allow for insertion site monitoring.



High MVTR

Provides environment with optimal moisture management



Secure Adhesion

Adhesive dressing enables firm fixation, minimising device movement or premature displacement



Waterproof

Protective layer prevents liquid from passing through covering




Suitable for:

- Protection and securement of PICC lines
- Protection and securement of Ports
- Protection and securement of PIVC lines
- Protection and securement of CVC lines

REF	Size	Type	Sterility	Pcs/Unit	Units/Ctn
ISB403	10cm x 12cm	Bordered	Sterile	100	12
ISB404	11cm x 15cm	Bordered	Sterile	100	12

Component details →

IV Film Dressing Range




ASB301



ASB302



ASB303



ASB304



ASB312



ASB402



ASB403/4



ASB202/3



IVD002

REORDER	SIZE	SUITABLE FOR	STERILITY	PCS/UNIT	UNITS/CTN
ASB301	5cm x 5.7cm	Paediatric, PIVC	Sterile	100	10
ASB302	6.5cm x 7cm	PIVC	Sterile	100	20
ASB303	7cm x 8.5cm	PIVC	Sterile	100	18
ASB304	7cm x 9cm	PIVC	Sterile	100	12
ASB312	8.5cm x 10.5cm	PIVC, CVC	Sterile	100	16
ASB402	6cm x 7cm	Paediatric, PIVC	Sterile	100	24
ASB403	10cm x 12cm	PICC, PORTS, PIVC, CVC	Sterile	100	12
ASB404	11cm x 15cm	PICC, PORTS, PIVC, CVC	Sterile	100	12
ASB202	6cm x 7cm	PICC, PORTS, PIVC, CVC	Sterile	50	10
ASB203	10cm x 12cm	PICC, PORTS, PIVC, CVC	Sterile	25	10
IVD002	6cm x 7cm	PICC	Sterile	50	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

