

## Non-Woven Roll

Asguard Non-Woven Tape Roll is a “cut to length” self-adhesive, non-woven fabric fixation tape. Highly comfortable, this tape is ideal for fixation for difficult contour areas without movement restriction. This versatile tape roll is ideal for window dressing as well as securing and aiding tubes, catheters, and other medical devices remain firmly in place, providing patients with the stability and comfort they need during their recovery process.



### Soft Non-Woven Material

Gentle on the skin, aiding to minimise the risk of irritation or discomfort as well as reduces edge lifting



### Cut to Length

Measuring grid backing allows for quick and accurate cutting to required size



### Flexible

Ideal for use on areas difficult to dress (e.g. heel)



### Suitable for:

- Fixation over joints and body contours
- Fixation of dressings and access devices
- Skin protective layer under more aggressive adhesives

REF	Size	Sterility	Pcs/Unit	Units/Ctn
ISD020	2.5cm x 10m	Non-Sterile	1	30
ISD021	5cm x 10m	Non-Sterile	1	24
ISD022	7.5cm x 10m	Non-Sterile	1	24
ISD023	10cm x 10m	Non-Sterile	1	12
ISD024	15cm x 10m	Non-Sterile	1	12
ISD025	20cm x 10m	Non-Sterile	1	6

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

