

# Asguard<sup>®</sup>

**SILICONE BORDER LITE**








# Asguard Silicone Border Lite


Asguard Silicone Border Lite is a 3 Layer, thin, ultra-flexible foam dressing featuring a gentle silicone adhesive. Designed for a wide variety of low-exuding wounds, it's ideal whenever a low-profile foam dressing is required. Combining effective absorption with enhanced mobility, Asguard Silicone Border Lite supports wound management without compromising comfort or flexibility.



## Why Choose Asguard Silicone Border Lite?

	<b>ATRAUMATIC REMOVAL</b>	Designed to reduce discomfort during dressing changes, the soft silicone adhesive minimizes pain and tissue damage.
	<b>SKIN FRIENDLY PERFORMANCE</b>	Ideal for fragile or sensitive skin, offering protection against trauma and making it a smart choice for elderly and paediatric patients.
	<b>EXTENDED SILICONE BORDER</b>	The increased size and width of the silicone border ensures longer wear time and reduces the frequency of dressing changes.
	<b>COMFORT IN MOTION</b>	Flexible and soft-edged silicone allows the dressing to conform to body contours and mobile areas such as knees and elbows, reducing the risk of rolling or detachment.
	<b>PROTECTIVE BARRIER</b>	Provides a barrier against external contaminants and protects the wound during daily activities such as showering.

## Multi-Layer Design



**Absorbent  
Foam Core**



**Water  
Resistant Film**



**Asguard  
Silicone**




# Secure, Gentle, Effective

## Asguard Silicone Border Lite

Compared to traditional adhesives or contact layers, Asguard Silicone high-adhesive technology offers a balance of secure fixation, gentle removal, and clinical versatility. It is especially beneficial for:

- Chronic wounds
- Fragile or sensitive skin
- Pediatric and geriatric patients
- Post-surgical sites



### Silicone Adhesive Layer

High Adhesion featuring a larger Silicone Border for improved securement and extended wear time



### Low-Profile Design

Thin, smooth layers reduce bulk and friction, ensuring compatibility under compression or clothing



### Absorbent Foam Core

Manages exudate efficiently



### Water & Shower Resistant Film

Designed to reduce dressing changes and maintain integrity during routine hygiene and mobility



### Available in Multiple Sizes

Square and rectangular formats for different wound types and locations



\* The specified goods should only be used under the supervision of a healthcare professional.

\* If you are concerned about your wound, consult a healthcare professional



REF	SIZE	PCS/UNIT	UNITS/CTN
<b>AWS20505</b>	5cm x 5cm	50	10
<b>AWS05125</b>	5cm x 12.5cm	50	10
<b>AWS7575</b>	7.5cm x 7.5cm	25	10
<b>AWS1010</b>	10cm x 10cm	100	20
<b>AWS1515</b>	15cm x 15cm	100	18