

## SurgiFlex Post-Operative Silicone Dressing

Surgiflex Si Silicone Dressing is a highly absorbent all-in-one island dressing that protects both skin and wound tissue. Absorbing excess exudate, this dressing provides the benefits of modern wound healing and encourages less frequent dressing changes. SurgiFlex Si is ideal for post-surgical procedures where exudate monitoring and management is key.



### Silicone Contact Layer

Gentle on skin and newly formed tissue



### Exudate Management

Multi-layer construction to increase exudate absorption



### Conforms to Body Contours

Spiral embossed pattern provides increased conformability



### Suitable for:

- Fragile, sensitive or reactive skin
- Protection of newly formed tissue
- Post-operative wounds and incisions

REF	Dressing Size	Padding Size	Sterility	Pcs/Unit	Units/Ctn
<b>SGF6070</b>	6cm x 7cm	3.5cm x 4.5cm	Sterile	10	10
<b>SGF100100</b>	10cm x 10cm	6.5cm x 6.5cm	Sterile	10	10
<b>SGF100200</b>	10cm x 20cm	5cm x 15cm	Sterile	10	10
<b>SGF100250</b>	10cm x 25cm	5cm x 20cm	Sterile	10	10
<b>SGF100300</b>	10cm x 30cm	5cm x 25cm	Sterile	10	10
<b>SGF100350</b>	10cm x 35cm	5cm x 30cm	Sterile	10	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

