

Asguard[®]

SURGIFLEX POST-OPERATIVE SILICONE DRESSING



SURGIFLEX POST-OPERATIVE SILICONE DRESSING


For over 30 years, Sentry Medical has been helping health care professionals provide exceptional care to their patients through a modern and quality driven portfolio.

The combination of silicone interface with superabsorbent fibres provides superior absorbency for newly forming tissue

Surgiflex Post-Operative Silicone Dressing is a highly absorbent all-in-one island dressing that protects both skin and wound tissue.

Absorbing excess exudate, this dressing provides the benefits of modern wound healing and minimises dressing changes. Surgiflex is ideal for post-surgical procedures where exudate monitoring and management is key.

Spiral pattern provides flexibility over joints



Silicone atraumatic interface provides gentle adherence

Shower proof backing allowing versatile patient activity

Did You Know?



Surgiflex Post-Operative Silicone Dressing is a combination of a silicone wound and skin contact layer, a superabsorbent polymer core, backed with a transparent polyurethane film outer layer



The combination of fluid absorption capacities with in-built moisture vapour transmission (MVTR) makes it suitable for the management of wounds with moderate to heavy levels of exudate



In the presence of moisture, the non-adherent property of the silicone layer of Surgiflex Post-Operative Silicone dressing prevents adherence to the wound, newly formed tissue and surrounding fragile, or sensitive skin



The outer backing of Surgiflex Post-Operative Silicone Dressings outer backing is water resistant, making activities of daily living, such as showering, easy to accommodate

Surgiflex Post-Operative Silicone Dressing Features:

ATRAUMATIC INTERFACE



The silicone layer of Surgiflex dressing prevents adherence to the wound and newly formed tissue and is gentle on surrounding fragile, or sensitive skin

WATER RESISTANCE



The outer backing of Surgiflex Post-Operative Silicone Dressing is water resistant, making activities of daily living, such as showering, easy to accommodate

FLEXIBLE & CONFORMABLE



The addition of a spiral embossed pattern provides conformability across the dressing for added comfort on joints and moving

EXUDATE MANAGEMENT



Multi-layer construction and fluid absorption capacities makes it manageable for all types of superficial wounds with light to moderate to heavy exudate production

TRANSPARENT DRESSING



Surgiflex Post-Operative Silicone Dressing is backed with a transparent polyurethane film outer layer which enables easy inspection of the peri-wound

VARIOUS SIZES



Can be used on a wide range of wounds and is available in multiple sizes



REF	DRESSING SIZE	PADDING SIZE	PCS/UNIT	UNITS/CTN
SGF6070	6cm x 7cm	3.5cm x 4.5cm	10	10
SGF100100	10cm x 10cm	6.5cm x 6.5cm	10	10
SGF100200	10cm x 20cm	5cm x 15cm	10	10
SGF100250	10cm x 25cm	5cm x 20cm	10	10
SGF100300	10cm x 30cm	5cm x 25cm	10	10
SGF100350	10cm x 35cm	5cm x 30cm	10	10