

Transparent Film Dressings

The Asguard Film Dressing is a transparent polyurethane adhesive film with moisture vapour permeability which adheres to the skin and forms a flexible and comfortable layer. Asguard Film Dressings conform to the skin and allows easy visibility to wound healing. They are ideal for the management of all types of insertion sites to provide a waterproof environment and/or to keep primary dressings in place.



Waterproof

Outer layer protects site from fluid penetration, allowing patient to shower with the dressing in-situ



Transparent Film

Allows visualisation of covered tissue



Easy Peel Frame

Assists with smooth and easy application



Suitable for:

- Wounds with little or no exudate
- Protection and securement of IV lines
- Dressing protection and reinforcement

REF	Size	Sterility	Pcs/Unit	Units/Ctn
ISB001	4cm x 5cm	Sterile	50	10
ISB002	6cm x 7cm	Sterile	50	10
ISB003	10cm x 12cm	Sterile	25	10
ISB004	15cm x 20cm	Sterile	25	10
ISB1025	10cm x 25cm	Sterile	20	10
ISB2030	20cm x 30 cm	Sterile	10	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

