

Transparent Film Island Dressing

Asguard Transparent Film Island Dressings provide a low adherent and breathable absorbent pad permeable to moisture. The waterproof film assists with unnecessary dressing changes. The adhesive is skin-friendly and provides reliable fixation, yet is easy to remove, minimising trauma to the healing tissue.



Flexible

Ideal for use on areas difficult to dress (e.g. heel)



Waterproof

Outer layer protects site from fluid penetration, allowing patient to shower with the dressing in-situ



Low Adherent Pad

Comfortable and absorbent with minimal pain on removal



Suitable for:

- Post-operative wounds
- Minor cuts and abrasions
- Wounds with a low level of exudate

REF	Dressing Size	Padding Size	Sterility	Pcs/Unit	Units/Ctn
ISC001	4cm x 5cm	2cm x 2.5cm	Sterile	50	10
ISC002	6cm x 7cm	3cm x 4cm	Sterile	50	10
ISC0610	6cm x 10cm	2.5cm x 6cm	Sterile	50	4
ISC0910	9cm x 10cm	4.5cm x 6cm	Sterile	25	10
ISC0915	9cm x 15cm	4.5cm x 10cm	Sterile	25	10
ISC003	10cm x 12cm	4.5cm x 6cm	Sterile	25	10
ISC028	10cm x 20cm	5cm x 15cm	Sterile	25	10
ISC025	10cm x 25cm	5cm x 20cm	Sterile	25	10
ISC026	10cm x 30cm	5cm x 25cm	Sterile	25	10
ISC0935	10cm x 35cm	5cm x 30cm	Sterile	30	8
ISC004	15cm x 20cm	8cm x 10.5cm	Sterile	25	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

