

## Undercast Padding

Asguard Undercast Paddings are made from soft and cushioning material, designed to provide comfort and protection. The padding is intended to be applied to and/or compress a local area of the body for various preventative and therapeutic applications such as in compression bandaging, padding under casts and splints and padding for anatomical depressions.



**Conforms to Body Contours**  
Allows flexibility without restriction of movement



**Soft Cushion**  
Provides comfortable cushion layer



**Air Permeable**  
Reduces the likelihood of sweating



- Suitable for:**
- Synthetic casts
  - Padding under casts, splints, or for anatomical depressions
  - Compression bandaging

| REF       | Size         | Type    | Sterility   | Pcs/Unit | Units/Ctn |
|-----------|--------------|---------|-------------|----------|-----------|
| 10105711* | 5cm x 2.7m   | Viscose | Non-Sterile | 12       | 30        |
| 10105712  | 7.5cm x 2.7m | Viscose | Non-Sterile | 12       | 20        |
| 10105713  | 10cm x 2.7m  | Viscose | Non-Sterile | 12       | 15        |
| 10105714  | 15cm x 2.7m  | Viscose | Non-Sterile | 6        | 20        |
| 10105731  | 5cm x 2.7m   | Cotton  | Non-Sterile | 12       | 30        |
| 10105732  | 7.5cm x 2.7m | Cotton  | Non-Sterile | 12       | 20        |
| 10105733  | 10cm x 2.7m  | Cotton  | Non-Sterile | 12       | 15        |
| 10105734  | 15cm x 2.7m  | Cotton  | Non-Sterile | 6        | 20        |

\*Indent item, minimum order quantities apply.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

