

# bodichek®

## Infrared Contactless Thermometer



### Features

CE certified

Features two modes – forehead and object

Super fast - one second response time

Temperature range 34.0°C~43.0°C(93.2°F~109.4°F)

Fever alarm when temperature is over 37.8°C (100.4°F)

Switches easily between Celsius and Fahrenheit readings

Memory function with real time clock

Stores up to ten memory readings each for forehead and objects

Jumbo digital display with back light

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
13011237	Bodichek® Infrared Contactless Thermometer	100	N/A	100



Aaxis Pacific  
7A Bessemer St. Blacktown  
NSW 2148 AUSTRALIA

T. +61 2 9881 3333  
F. +61 2 9881 3322  
E. info@aaxispacific.com.au

aaxis.com.au  
V3.0