



# Bodichek® Digital Water Proof Thermometers



13011226



13011229



## Features

- 
- Waterproof IPx7
- 
- Complies with ACCC and TGA button battery safety guidelines
- 
- Measurement range: 32° C – 41.9° C
- 
- Large, clear, easy to read digital display
- 
- Automatically turns off after 10 minutes of inactivity
- 
- Convenient beeper signal indicating thermometer can be removed from reading
- 
- Last temperature memory recalled for later reference if needed
- 



13014001

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
13011226	Bodichek® Digital Thermometer, Water Proof, Rigid Tip	12	12	144
13011229	Bodichek® Digital Thermometer, Water Proof, Flexible Tip	10	10	120
13014001	Bodichek® Digital Thermometer Probe Cover 100's	12	12	120