



Bodifast® Tubular Retention Bandage

Bodifast® is a high quality seamless white tubular bandage which securely holds wound dressings in place without restricting movement or creating a tourniquet effect. Due to their versatility, Bodifast® bandages are primarily used to securely hold wound dressings in place, as a light dermatological covering for eczema sufferers (wet wrap) and as an under cast/orthopaedic device netting.



Suggested applications

- Securely holds wound dressings in place
- Light dermatological covering for eczema (wet wrap)
- Under cast/orthopaedic device netting
- Bandage replacement solution for most wound types
- Skin covering for any part of the body

Features

- Seamless bandage design (no pins or tapes required)
- Light elasticity prevents compression
- Stretchable, allowing freedom of movement
- Breathable, preventing wound maceration
- Colour coding design to easily distinguish between sizes
- Hypoallergenic (latex free)

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
12823510	Bodifast® Tubular Retention Bandage 3.5cm x 10m Red	20	20	80
12823511	Bodifast® Tubular Retention Bandage 5cm x 10m Green	16	16	64
12823512	Bodifast® Tubular Retention Bandage 7.5cm x10m Blue	12	12	48
12823513	Bodifast® Tubular Retention Bandage 10cm x 10m Yellow	8	8	32
12823515	Bodifast® Tubular Retention Bandage 20cm x 10m Purple	4	4	16
12823500	Bodifast® Tubular Retention Bandage 3.5cm x 2m Red	12	12	72
12823501	Bodifast® Tubular Retention Bandage 5cm x 2m Green	12	12	72
12823502	Bodifast® Tubular Retention Bandage 7.5cm x 2m Blue	12	12	72
12823503	Bodifast® Tubular Retention Bandage 10cm x 2m Yellow	12	12	72
12823505	Bodifast® Tubular Retention Bandage 20cm x 1m Purple*	12	12	72

*Purple size comes in 1 metre