



Bodigrip® Multi-Purpose Tubular Support Bandage

Bodigrip® tubular compression bandages have been specially designed to provide effective support for sprains, strains, weak joints, soft tissue injuries, most types of oedema and post burn scarring. Bodigrip® tubular compression bandages are Type 3A & 3B – light to moderate compression. The bandage distributes pressure evenly and fits to the contours of the body, allowing freedom of movement and improving comfort. Bodigrip® bandages are made with elastane with excellent stretch and recovery properties and encased in two layers of cotton to protect sensitive and fragile skin.



Features

Available in a range of sizes to suit all limbs and joints

Suitable for pressure dressings

Easy to apply with no need for pins or tapes

Provides continuous support with low level pressure



Bodigrip® Multi-Purpose Tubular Support Bandage

User Instructions:

Measure around the widest part of the affected area and select the appropriate size

Cut bandage to twice the required length for the affected area of the limb

Pull over the affected area

Double back over limb

Pull top layer about 3cm past bottom layer



PRODUCT CODE (Natural)	PRODUCT CODE (Flesh)	SIZE	APPLICATION	BANDAGE WIDTH	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
13023100	13023120	A	Infant feet/arms	4.5cmx10m	1	N/A	40
13023101	13023121	B	Small hands/arms	6.5cmx10m	1	N/A	28
13023102	13023122	C	Medium arms, small ankles	6.75cmx10m	1	N/A	28
13023103	13023123	D	Large arms, medium ankles, small knees	7.5cmx10m	1	N/A	24
13023104	13023124	E	Large ankles, medium knees, small thighs	8.75cmx10m	1	N/A	20
13023105	13023125	F	Large knees, medium thighs	10cmx10m	1	N/A	20
13023106	13023126	G	Large thighs	12cmx10m	1	N/A	16
13023107	13023127	J	Small trunks	7.5cmx10m	1	N/A	12
13023108	13023128	K	Medium trunks	21.5cmx10m	1	N/A	8
13023109	13023129	KL	Large trunks	27cmx10m	1	N/A	8
13023110	13023130	L	X-Large trunks	32.5cmx10m	1	N/A	6
13023111	13023131	M	Extra X-Large trunks	36.5cmx10m	1	N/A	6