



Basic Isolation & X-Ray Gowns

S+M® Isolation Gowns are a single non-woven polypropylene disposable gown equally suitable for patients or staff with light fluid protection. The gowns have a large fit size and come in either full body and arm coverage or sleeve version.



Features

Non-woven lightweight polypropylene (PP)

Low lint materials and fully recyclable

Neck and waist ties, generously sized for full coverage and ease of movement

30 gsm for Isolation Gowns with Sleeve and Cuffed Sleeve, 14 gsm for Isolation Gown Apron – polyethylene (PE)

Available in various colours for ward purpose identification

PRODUCT CODE	PRODUCT DESCRIPTION	TOTAL LENGTH (mm)	SHOULDER WIDTH (mm)	SEE THROUGH	SKU/TU (MOQ)	SKU/Inner	SKU/CTN
10125311	S+M® Gown Isolation X-Ray M Blue	1050	740	N	100	25	100
10125313	S+M® Gown Isolation X-Ray L Blue	1300	800	N	100	25	100
10125312	S+M® Gown Isolation Sleeve Yellow	1050	710	Y	50	10	50
10125314	S+M® Gown Isolation Sleeve Light Blue	1050	710	Y	50	10	50
20126112	S+M® Disposable Apron PE 81x135cm	1350	810	Y	1000	100	1000
PA61001Y	Suresafe® Basic Isolation Gown, Yellow	1150	1370	Y	100	N/A	100