


# Bed Roll & Stretcher Covers



## Suggested Applications

Ideal for use in hospitals, medical centres, GP clinics, physiotherapists, aged care and radiology units

Suitable for: bench/stretcher/table sheet; head/pillow cover; underlay; tray liner

S+M® Bed Roll completely covers the majority of NSW Ambulance stretchers

## Features

Strong, absorbent and water-resistant

Perforated at 50 cm intervals, no scissors required to cut to length

Dual layer draws and locks moisture away from patient

Made from 24gsm paper with a polyethylene backing

72 cm width covers the entire width of examination couch and meets the new guidelines of RACGP exam couches - wider than competitors standard of 53.5cm

Easily disposable

Cost effective

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10125517	S+M® Bed Cover Bed Roll 72cmx80m Light Blue	6	N/A	6
10125518	S+M® Bed Cover Bed Roll 72cmx80m White	6	N/A	6
10130220	S+M® Bed Cover Roll Dispenser	1	N/A	1
10125513	S+M® Bed Cover Stretcher Cover 70x240cm Navy Blue	100	10	100
10125514	S+M® Bed Cover Stretcher Cover 70x240cm White	100	10	100
10025513	AAXIS® Bed Stretcher Cover 70x240cm Navy Blue	120	10	120
10025514	AAXIS® Bed Stretcher Cover 70x240cm White	120	10	120