

## CHG Swab Stick

CHG Swab Sticks consist of a non-woven pledget adhered to a plastic stick, used for cleaning the skin before a procedure. The pledget is saturated with 2% chlorhexidine gluconate and 70% isopropyl alcohol solution.



**Durable**  
Built with a sturdy but flexible plastic stick



**Individually Wrapped**  
Individually sealed in easy tear foil sachets



**2% CHG, 70% IPA**  
Saturated with a 2% chlorhexidine gluconate and 70% isopropyl alcohol solution



**Suitable for:**

- Hospital and medical
- Topical cleansing prior to procedure

REF	Size	Quantity	Sterility	Pcs/Unit	Units/Ctn
10012610	13 cm	50 pcs	Non-Sterile	50	20

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



**Sentry Medical Pty Ltd**  
1300 995 999 | [sales@sentrymedical.com.au](mailto:sales@sentrymedical.com.au)  
[www.sentrymedical.com.au](http://www.sentrymedical.com.au)

