

CHG Swab Stick

CHG Swab Sticks consist of a non-woven pledget adhered to a plastic stick, used for cleaning the skin before a procedure. The pledget is saturated with 2% chlorhexidine and 70% isopropyl alcohol solution.



Durable

Built with a sturdy but flexible plastic stick



Individually Wrapped

Individually sealed in easy tear foil sachets



70%

70% Isopropyl Alcohol

Saturated with a 2% chlorhexidine and 70% isopropyl alcohol solution



Suitable for:

- Hospital and medical
- Topical cleansing prior to procedure

REF	Size	Quantity	Sterility	Pcs/Unit	Units/Ctn
10012610	13.7 cm	50 pcs	Non-Sterile	50	20

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.