

CHG Swab Sticks and Swabs

Bodichek CHG Swab Sticks and Swabs consist of a non-woven swab saturated with chlorhexidine gluconate and isopropyl alcohol solution, designed to be used for cleaning the skin before a procedure (alcohol-free option available). The Bodichek Swab Stick pledgets are adhered to a sturdy but flexible plastic stick for ease of application.



Ideal for Skin Prep

For cleaning localised body surface prior to a procedure or device application




Easy Tear Sachets

Sealed in easy-to-tear foil sachets for user convenience and ease of packaging removal



Swab Sticks offer a Flexible and Durable Stick with Isolated Tip

The pledgets ensure clinicians can apply solution without making physical contact with the skin



10012610

Suitable for:

- Hospital and medical
- Topical cleansing prior to procedure

REF	Description	Solution	Size	Sterility	Pcs/Unit	Units/Ctn
10012610	Swab Stick	2% CHG, 70% IPA	13cm	Non-Sterile	50	20
10012611*	Swab Stick	2% CHG, 70% IPA	14.8cm	Non-Sterile	50	20
10012612*	Swab Stick	0.5% CHG, 70% IPA - Red Tint	14.8cm	Non-Sterile	50	20
10012613*	Swab Stick	2% CHG, Alcohol-Free	14.8cm	Non-Sterile	50	20
10012614*	Swab Stick	2% CHG, 70% IPA	7.4cm	Non-Sterile	60	20
10012615*	Swab Stick x 3	2% CHG, 70% IPA	7.4cm	Non-Sterile	60	20
10012620*	Swab	2% CHG, 70% IPA	4.5cm x 8.5cm	Non-Sterile	100	100
10012621*	Swab	2% CHG, 70% IPA	6cm x 6cm	Non-Sterile	200	50

*Indent item, minimum order quantities apply

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

