

Bodichek®

Digital Thermometer Probe Covers

Bodichek Digital Thermometer Probe Covers are single-use probe covers designed to cover a thermometer probe when in use, acting as a physical barrier to help prevent cross-contamination between patients.



Cross-Contamination Protection

Acts as a barrier to prevent cross-contamination and the proliferation of germs




Disposable

Single-use probe covers to help reduce the risk of spreading infections between patients



Not Made with Natural Rubber Latex

For users with allergies or sensitivities



Perforated dispenser box for ease of use.

Suitable for: - Hospital and medical

REF	Sterility	Pcs/Unit	Unit/Inner	Units/Ctn
13014001	Non-sterile	100	12	120

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

