

Elastic Adhesive Bandages

Bodichek Elastic Adhesive Bandages are flexible cotton bandages with zinc oxide adhesive on one side, allowing them to provide compression and long-lasting support for injured or strained muscles, joints, or ligaments.



Conforms to Body Contours

Allows flexibility without restriction of movement



High Shear Adhesive

Secure hold, ensuring that the bandage stays in place during treatment period



Alignment Reference Line

Central yellow line assists in ensuring proper placement during the application process



- Suitable for:**
- Support for soft tissue injuries such as strains and sprains
 - Assisting with injury prevention
 - Compression to alleviate swelling

REF	Size	Feature	Sterility	Pcs/Unit	Units/Ctn
10105321	5cm x 1.8m	Unstretched	Non-Sterile	12	24
10105322	7.5cm x 1.8m	Unstretched	Non-Sterile	12	16
10105323*	10cm x 1.8m	Unstretched	Non-Sterile	12	12

*Indent item, minimum order quantities apply

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

