

## Fabric Strips

Bodichek Fabric Strips are intended for the protection and management of minor superficial wounds. The strip helps protect against water, dirt, and everyday contamination while remaining comfortable on the skin. The strip may also be layered to provide additional coverage and light compression where clinically appropriate.



### Breathable

Allows airflow to support skin integrity and patient comfort during wear




### Protective Barrier

The strip assists in shielding the wound site from environmental contaminants and particulates



### Gentle on Skin

Low-adherent absorbent pad minimises trauma to the wound bed and supports atraumatic removal



13010001



- Suitable for:**
- Minor lacerations, cuts and abrasions
  - High-mobility anatomical areas (e.g. fingers, elbows, knees)
  - Securing primary wound dressings

| REF      | Description   | Size        | Sterility | Pcs/Box | Bxs/Shipper | Shippers/Ctn |
|----------|---------------|-------------|-----------|---------|-------------|--------------|
| 13010001 | Fabric Strips | 20mm x 72mm | Sterile   | 25      | 12          | 144          |
| 13010002 | Fabric Strips | 20mm x 72mm | Sterile   | 50      | 8           | 96           |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

