

Plastic Strips

Bodichek Plastic Strips provide reliable protection for minor wounds and abrasions. Each strip features a low-adherent absorbent pad, and a flexible plastic backing with breathable perforations. The strip helps protect against water, dirt, and everyday contamination while remaining comfortable on the skin.



Breathable

Allows airflow to support skin integrity and patient comfort during wear



Protective Barrier

The strip assists in shielding the wound site from environmental contaminants and particulates



Gentle on Skin

Low-adherent absorbent pad minimises trauma to the wound bed and supports atraumatic removal



13010003

- Suitable for:**
- Minor lacerations, cuts and abrasions
 - Superficial wounds

REF	Size	Sterility	Pcs/Box	Bxs/Shipper	Shippers/Ctn
13010003	20mm x 72mm	Sterile	50	8	96
13010004	20mm x 72mm	Sterile	100	8	80

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

