

Flat Rolls Tyvek®

Produced with Tyvek® on one side and a PET/PP blend film on the other, Bolsaplast Flat Rolls offer extremely low fibre shedding and are resistant to film tearing. The rolls are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.



**Colour Change Indicator**  
For plasma sterilisation (H2O2) printed under seal



**Non-Toxic Ink**  
No damage to instruments



**Medical Grade Construction**  
Offers excellent peelability and is resistant to film tearing



**Compliance**  
Conforms to EN868-2/3/4/5 and ISO11607-1



**Suitable for:** - Sterilising devices and equipment

REF	Size	Sterility	Pcs/Unit	Units/Ctn
MRT1171	7.5cm x 70m	Non-sterile	1	1
MRT1172	10cm x 70m	Non-sterile	1	1
MRT1174	15cm x 70m	Non-sterile	1	1
MRT1175	20cm x 70m	Non-sterile	1	1
MRT1176	25cm x 70m	Non-sterile	1	1
MRT1177<<	30cm x 70m	Non-sterile	1	1
MRT1178<<	40cm x 70m	Non-sterile	1	1
MRT1179<<	50cm x 70m	Non-sterile	1	1

<< Indent Item

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

