

Heat Sealable Pouches

Produced with medical grade paper on one side and a PET/PP blend film on the other, Bolsaplast Heat Sealable Pouches offer extremely low fibre shedding and are resistant to film tearing. The pouches are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.



Colour Change Indicator
Chemical indicators for Steam/EO/formaldehyde



Non-Toxic Ink
No damage to instruments



Medical Grade Construction
Offers excellent peelability and is resistant to film tearing



Suitable for: - Sterilising devices and equipment

Reorder	Size	Pcs/Unit	Units/Ctn
MBO1025	7.5cm x 15cm	2700	1
MBO1000	7.5cm x 21cm	2500	1
MBO1001	7.5cm x 28cm	2400	1
MBO1003	10cm x 20cm	2400	1
MBO1005	10cm x 28cm	2400	1
MBO1006	10cm x 42cm	2400	1
MBO1022<<	10cm x 58cm	2000	1
MBO1007	15cm x 21cm	2000	1
MBO1008	15cm x 28cm	2000	1
MBO1009	15cm x 32cm	2000	1
MBO1010	15cm x 36cm	2000	1

Reorder	Size	Pcs/Unit	Units/Ctn
MBO1021	15cm x 40cm	1200	1
MBO1011	20cm x 32cm	1800	1
MBO1012<<	20cm x 36cm	1800	1
MBO1013	20cm x 42cm	1200	1
MBO1023<<	20cm x 50cm	1200	1
MBO1014	25cm x 36cm	1200	1
MBO1015	25cm x 42cm	1200	1
MBO1016	25cm x 50cm	1200	1
MBO1017	30cm x 42cm	600	1
MBO1024<<	30cm x 58cm	700	1
MBO1019<<	40cm x 58cm	600	1

<< Indent Item

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.