

Sterilisation Flat Reel

Produced with medical grade paper on one side and a PET/PP blend film on the other, Bolsaplast Sterilisation Flat Reels offer extremely low fibre shedding and are resistant to film tearing. The reels are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.



Colour Change Indicator
Chemical indicators for Steam/EO/formaldehyde



Non-Toxic Ink
No damage to instruments



Medical Grade Construction
Offers excellent peelability and is resistant to film tearing



Compliance
Conforms to EN868-2/3/4/5 and ISO11607-1



Suitable for: - Sterilising devices and equipment

REF	Size	Sterility	Pcs/Unit	Units/Ctn
MRP1040	5cm x 200m	Non-sterile	1	1
MRP1041	7.5cm x 200m	Non-sterile	1	1
MRP1042	10cm x 200m	Non-sterile	1	1
MRP1043 <<	12.5cm x 200m	Non-sterile	1	1
MRP1044	15cm x 200m	Non-sterile	1	1
MRP1045	20cm x 200m	Non-sterile	1	1
MRP1046	25cm x 200m	Non-sterile	1	1
MRP1047	30cm x 200m	Non-sterile	1	1
MRP1048	40cm x 200m	Non-sterile	1	1

<< Indent Item

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.