

Sterilisation Gusset Reels

Produced with medical grade paper on one side and a PET/PP blend film on the other, Bolsaplast Sterilisation Gusset Reels offer extremely low fibre shedding and are resistant to film tearing. The reels are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.



Colour Change Indicator
Chemical indicators for Steam/EO/formaldehyde



Non-Toxic Ink
No damage to instruments



Medical Grade Construction
Offers excellent peelability and is resistant to film tearing



Compliance
Conforms to EN868-2/3/4/5 and ISO11607-1



Suitable for: - Sterilising devices and equipment

REF	Size	Sterility	Pcs/Unit	Units/Ctn
MRF1060	10cm x 4cm x 100m	Non-sterile	1	1
MRF1061	15cm x 4cm x 100m	Non-sterile	1	1
MRF1062	20cm x 6cm x 100m	Non-sterile	1	1
MRF1063	25cm x 6cm x 100m	Non-sterile	1	1
MRF1064<<	30cm x 8cm x 100m	Non-sterile	1	1
MRF1065<<	40cm x 8cm x 100m	Non-sterile	1	1

<< Indent Item

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

