



Strong and Durable



Various Sizes



Transparent

Dust Cover Bags provide an ideal storage environment for instruments and trays that have undergone the sterilisation process, and need protection from dust whilst in storage or transportation. These bags can be heat sealed using an impulse bar sealer.

HEAT SEALABLE

Provides fast and easy closing with impulse bar sealer

TRANSPARENT

Clear for ease of viewing contents

HIGH PROTECTION

Maintains integrity of sterile equipment during storage



Suitable for:

Storage of sterilised instruments and trays

Extends sterile tenability

Reorder	Size	Colour	Pcs/Unit	Units/Ctn
DC001	15.2cm x 25.4cm	Clear	100	20
DC002	20.3cm x 30.5cm	Clear	100	15
DC003	25.4cm x 38.1cm	Clear	100	15
DC004	30.5cm x 50.8cm	Clear	100	15
DC005	40.6cm x 50.8cm	Clear	100	5
DC006	45cm x 76.2cm	Clear	100	4
DC007	61cm x 76.2cm	Clear	100	4