



AAMI
Level 2



Impervious



Various
Colours

Owear Impervious Gowns provide effective protection against fluids with the front and back of the gown being impervious. Isolation Gowns comply with ISO6529 standards, as well as meet AAMI liquid barrier performance levels.

IMPERVIOUS

An impervious layer prevents fluid to pass through

VELCRO NECK AND TIE WAIST


Velcro neckline and back ties to easily don and doff

ELASTIC CUFF


Durable and easy to don and doff



Blue



Yellow



Purple
(Oncology)

Suitable for:

Hospital and medical

Laboratory

Reorder	Size	Feature	Colour	Pcs/Unit	Units/Ctn
GN004P	Regular	AAMI Level 2	Purple (Oncology)	100	1
GN005P	Regular	AAMI Level 2	Blue	50	1
GN006P	Regular	AAMI Level 2	Yellow	50	1