

Sentry Medical Gauze Swab - 8 Ply



Low Linting



8ply Thickness



Various Sizes

100% cotton and are French folded to reduce fraying and linting. The gauze swabs are sterile and 8ply in thickness and ideal to be used as a wound dressing or absorbing fluid.

NON-FRAYING

French folded to prevent fraying

HIGHLY ABSORBENT

100% cotton providing high absorbency

Suitable for:

Cleaning wounds and soaking up fluid

Wound protection barrier

Sensitive skin



Reorder	Pcs	Size	Sterility	Colour	Pcs/Unit	Units/Ctn
GSR2223	3pcs	5cm x 5cm	Sterile	White	50	30
GSR2225	5pcs	5cm x 5cm	Sterile	White	50	24
GSR2322	2pcs	7.5cm x 7.5cm	Sterile	White	50	24
GSR2323	3pcs	7.5cm x 7.5cm	Sterile	White	50	24
GSR2325	5pcs	7.5cm x 7.5cm	Sterile	White	50	24
GSG2425	5pcs	10cm x 10cm	Sterile	Green	50	9
GSR2422	2pcs	10cm x 10cm	Sterile	White	50	16
GSR2423	3pcs	10cm x 10cm	Sterile	White	50	8
GSR2425	5pcs	10cm x 10cm	Sterile	White	50	8

