

ISB403 & ISB404 | AsGuard® Exceed I.V. Reinforced High MVTR Film Dressing

AsGUARD®
eXCEED I.V.



Water proof



High MVTR



Secure

The range of AsGuard® Exceed I.V. dressings are constructed from a High MVTR film, which allows for excess moisture to evaporate, combined with a laminated non-woven border for securement and enhanced fixation.

'HIGH TACK, SLOW RELEASE' ADHESION

Adhesive technology provides firm fixation, minimising device movement or premature displacement

EASY PEEL FRAME


Easy to grab tabs for simple peel back

LARGE TRANSPARENT WINDOW

Allows for optimal insertion site monitoring

HIGH MVTR

Provides environment with optimal moisture management



Suitable for:


Protection and securement of PICC lines


Protection and securement of Ports

Protection and securement of PIVC lines

Protection and securement of CVC line

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
ISB403	10cm x 12cm	Sterile	100	12
ISB404	11cm x 15cm	Sterile	100	12






ASB301



ASB302



ASB303



ASB304



ASB312



ASB402



ASB403/4



ASB202/3

REORDER	SIZE	SUITABLE FOR	STERILITY	PCS/UNIT	UNITS/CTN
ASB301	5cm x 5.7cm	Paediatric, PIVC	Sterile	100	10
ASB302	6.5cm x 7cm	PIVC	Sterile	100	20
ASB303	7cm x 8.5cm	PIVC	Sterile	100	18
ASB304	7cm x 9cm	PIVC	Sterile	100	12
ASB312	8.5cm x 10.5cm	PIVC, CVC	Sterile	100	16
ASB402	6cm x 7cm	Paediatric, PIVC	Sterile	100	24
ASB403	10cm x 12cm	PICC, PORTS, PIVC, CVC	Sterile	100	12
ASB404	11cm x 15cm	PICC, PORTS, PIVC, CVC	Sterile	100	12
ASB202	6cm x 7cm	PICC, PORTS, PIVC, CVC	Sterile	50	10
ASB203	10cm x 12cm	PICC, PORTS, PIVC, CVC	Sterile	25	10