

NoPress® Eye Pressure Shield



Flexible
Design



Comfortable
for Patient



Transparent
Shield

NoPress is designed to protect the anaesthetised patient's eyes from accidental pressure related injury due to positioning, draping and operational techniques or procedures.

EVEN PRESSURE

Single shield distributes even pressure across orbital ridges

FLEXIBLE DESIGN

Patented design allows flexation around nose and conformity to patient's facial shape

FOAM BORDER

No sharp down-facing edges which could injure underlying tissues if pressure is applied



Suitable for:

Hospital and medical

Operating theatre

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
IVG-03	One Size	Non-Sterile	20	5